

Believe in Innovation



Roller 20MC  
ESR in 20 Seconds





## FULLY AUTOMATED ANALYZER FOR THE DETERMINATION OF THE ESR WITH EXTERNAL NEEDLE FOR PEDIATRIC SAMPLES AND UNCAPPED TUBES

- Results in 20 seconds related to red cells aggregation
- Only 30 µl EDTA blood sample
- No reagents required
- Results expressed in mm/h
- High correlation with the Westergren method
- No influence of low hematocrit levels
- Use of the same CBC tubes
- Latex Calibration & Controls
- Smart cards
- Thermostated at 37°C
- LCD touch screen
- User-friendly software
- Connection to LIS
- Thermal printer
- External barcode reader

CE MARKED

### New generation CPS



Code SI R20-MC



### LATEX CONTROLS

### Smart CARDS

Roller 20MC works with the “universal” smart card, the same used for TEST1, Roller 20LC, Roller 20PN.

New Roller 20MC multicode card of 5.000 tests with possibility to activate batch of 500 tests at a time with a code number.

