



BOD Incubator:

JK Biomed is the most trusted Brand for BOD Incubator's offering various capacities with economical cost effective models, aesthetic appearance, wide optional features, Efficiently designed, Long Lasting Performance to meet daily routine requirement of the user under controlled temperature condition for determination of Bio chemical oxygen demand tests, Storage of media, Microbiology Studies, Plasma Storage, culture analysis studies, vaccine storage, germination process, food Industry, Biotech / Microorganism / Plant Life /fermentation studies, Bacterial Growth, Serum incubation & many more..

JK Biomed offers 8 different sizes of quality, durable & reliable BOD Incubators.



Salient Features:

- ✦ Exterior Powder coated or Stainless Steel model with Stainless Steel as interior..
- ✦ Double wall construction with insulated Body & door.
- ✦ Castor Wheel.
- ✦ Locking Facility.
- ✦ Dual Door with interior transparent.
- ✦ Air circulation with excellent uniformity & good accuracy levels obtained.
- ✦ Completely digital controlled models.
- ✦ Alarm Facility.
- ✦ Additional Safety Standby features.
- ✦ Durable hermetically sealed CFC free refrigeration system with long life & extended guarantee models.
- ✦ Durable & requires less maintenance.
- ✦ Reliable & easy to operate with long lasting performance.
- ✦ Different optional features available to suit user's requirement.
- ✦ Specialized & Custom built models also available to suit specific requirement of the user. For further details kindly contact us.



Sizes:

Models	Capacity	Interior Sizes	Exterior Sizes
	Ltrs	W X D X H	W X D X H
ES-BOD/100	100	18" x 16" x 22"	24" x 26" x 46"
ES-BOD/170	170	20" x 20" x 26"	26" x 30" x 52"
ES-BOD/275	275	24" x 22" x 32"	30" x 33" x 58"
ES-BOD/340	340	24" x 24" x 36"	30" x 35" x 64"
ES-BOD/450	450	28" x 26" x 38"	34" x 37" x 66"
ES-BOD/620	620	30" x 30" x 42"	36" x 41" x 72"
ES-BOD/800	800	32" x 32" x 48"	38" x 44" x 78"
ES-BOD/1000	1000	32" x 32" x 60"	38" x 44" x 90"